



► **Headquarter: Shenzhen China**

- ③ 3 R&D Centers
- ③ 3 Manufacturing Centers
- 70,000m² Production Area

- ③ Serving 16 Countries & Regions
- Serving 20,000+ Medical Institutions
- Serving 300 Million+ People

Member of CIRC, GIMC, 2020. December V1.0 Copyright © 2020 HEADWAY



Shenzhen Zhonghe Headway Bio-Sci & Tech Co., Ltd.

Add: 6/f, Building r1-b, South district, High-tech industrial park,
Shennan road, Shenzhen, 518057, China
Tel: +86-755-26619490 Fax: +86-755-26617554
Email: customer@headwaychina.com

www.headwaychina.com



HCBT-01 & ¹³C-UBT Kits



Diagnosis of H.Pylori infection





Founded in 1996 in Shenzhen, China, Headway is currently the leading manufacturer & supplier of H.Pylori diagnostic—¹³C & ¹⁴C urea breath test—systems in the world. We are committed to making it easier to pursue good health by leveraging our more than 20 years of innovation in breath tests.

Vision


To be the leader in global breath detection industry

Mission

Make the pursuit of health easy

Values

Innovation, Practicality, Responsibility, Cooperation



Barry Marshall and Robin Warren won 2005 Nobel Prize in medicine for first separating H.Pylori from a gastric mucosal specimen.

H.Pylori colonizes in the gastrointestinal track of nearly 50% of the world's population.

H.Pylori was classified as a class I human carcinogen by the World Health Organization in 1994.

H.Pylori represents a key factor in the aetiology of various

gastrointestinal diseases: Dyspepsia, Duodenal Ulcer, Gastric Ulcer, Chronic Gastritis, Gastric MALT lymphoma, MALT, Gastric Cancer.

H.Pylori is highly associated to the development of gastric cancer.

Various Testing Methods



UBT

Non-invasive, is perfectly suited for primary diagnosis and for post treatment follow-up of H.Pylori infection.



Histology

Invasive, several factors influence the accuracy of histology, such as site, size and number of biopsies, PPIs, antibiotics.



Culture

Invasive, less sensitive method. The cultivation in vitro requires particular transport medium, growth medium and incubation environment.



RUT

Invasive, several factors influence the accuracy, such as the density of bacteria present in the biopsy specimen, PPIs, bismuth, antibiotics and presence of blood.



Stool Antigen

Distasteful process of stool collection. Several factors influence the accuracy, like antibiotic, PPIs, bowel movement and upper gastrointestinal bleeding.



Serology

Low accuracy and the positive result may reflect previous rather than current infection, not useful after treatment.

Global guideline and consensus for UBT H.Pylori diagnosis



“UBT is the most reliable non-endoscopic test to document eradication of H. Pylori infection.”

American College of Gastroenterology



“The diagnostic accuracy of the UBT is >95% in studies. The UBT is an accurate, practical, and readily available test.”

European Helicobacter Study Group



“A distinct advantage of [UBT] is that it can also be used to determine treatment efficacy.”

Mayo Clinic Laboratories



UBT have been recommended as the most accurate noninvasive tests for diagnosis of H. Pylori infection and for confirmation of eradication after therapy.

Quest Diagnostics Laboratories



HCBT-01 features

Operation procedure of ¹³C-urea breath test

Headway provides a simple and clean solution for you and your family to protect stomach.

During the test, first you may simply breathe into baseline breath bag. After taking the Headway ¹³C-urea capsule and waiting for 30 minutes, you breathe into sample breath bag. If there is an infection, it can be detected in your exhaled breath with HCBT-01.



Collect baseline breath



Take one ¹³C-Urea capsule



Wait for 30 minutes



Collect sample breath



Use HCBT-01 to test the two bags of breath samples



Principle of HCBT-01

CO₂ has several characteristic absorption peaks in infrared band and there is a slight difference between the absorption peak positions of ¹³CO₂ and ¹²CO₂. By utilizing such slight difference, the HCBT-01 breath test analyzer is able to measure the concentrations of ¹³CO₂ and ¹²CO₂ respectively and then calculate the variation DOB (‰) of the abundance of ¹³C isotopes in the sample and baseline over the natural abundance of ¹³C, so as to determine whether the patient has been infected with H.Pylori.



Specifications

Operating Condition

Indoor (protected from direct sunlight, radiation)
Temperature: 15°C-30°C, avoiding sharp change in temperature
Elevation: ≤2000m
Relative humidity: ≤75%
Atmospheric pressure: 75kpa-106kpa

Storage and transportation

The storage environment should be protected against influences of temperature, humidity, dust, salty, corrosive gases and water
The instrument should be prevented from being inclined, inverted, shaken or drastically vibrated in transportation
Ambient temperature: -10°C- +60°C
Relative humidity: ≤ 75%

Physical Specifications

Dimensions
· Length: 540mm
· Width: 399mm
· Height: 333mm
Weight: 17Kg

Power supply

Power supply: AC110-240V
Maximum rated power consumption <250VA

